

# HEMO-ACCESS

## Single, Double and Triple Lumen Dialysis Catheter Sets

*Biometrix catheters are designed for temporary use for dialysis and hemofiltration in subclavian, jugular and femoral veins.*

Biometrix Dialysis Catheters Feature:

- Biocompatible thermosensitive catheter material softens after insertion for increased patient safety and comfort
- Color coded clamps and priming volume indication for easier catheter maintenance
- Rotatable suture wing provides secure yet flexible anchoring
- French size and length printed on hub for easier identification

### Catheters for Subclavian and Femoral Approach

3

#### Lumen Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length	Y-Introducer Needle	Quick Insertion Syringe
EG-7613	12Fr x 15 cm	60 cm	-	-
*EG-Y7613	12Fr x 15 cm	70 cm	-	-
EG-7713	12Fr x 20 cm	60 cm	-	-
*EG-Y7713	12Fr x 20 cm	70 cm	-	-
EG-I7613	12Fr x 15 cm	60 cm	+	-
EG-I7713	12Fr x 20 cm	60 cm	+	-
*EG-R7613	12Fr x 15 cm	60 cm	-	+
*EG-R7713	12Fr x 20 cm	60 cm	-	+

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.035" Marked Guidewire J-Flex Tip
- > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe
- > 2 Heparin Locks

(\*) EG-Y7613, EG-Y7713, EG-R7613, EG-R7713 contains a 0.035" Nitinol Guidewire

2

#### Lumen Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length	Y-Introducer Needle	Quick Insertion Syringe
EF-7613	12Fr x 15 cm	60 cm	-	-
*EF-Y7613	12Fr x 15 cm	70 cm	-	-
EF-7713	12Fr x 20 cm	60 cm	-	-
*EF-Y7713	12Fr x 20 cm	70 cm	-	-
EF-7813	12Fr x 24 cm	70 cm	-	-
*EF-Y7813	12Fr x 24 cm	70 cm	-	-
EF-I7613	12Fr x 15 cm	60 cm	+	-
EF-I7713	12Fr x 20 cm	60 cm	+	-
EF-I7813	12Fr x 24 cm	70 cm	+	-
*EF-R7713	12Fr x 20 cm	60 cm	-	+
*EF-R7813	12Fr x 24 cm	70 cm	-	+
*EF-9613	11Fr x 15 cm Rigid Tip	60 cm	-	-
*EF-9713	11Fr x 20 cm Rigid Tip	60 cm	-	-


Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.038" Marked Guidewire J-Flex Tip
- > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe
- > 2 Heparin Locks

(\*) EF-Y7613, EF-Y7713, EF-Y7813, EF-R7713, EF-R7813, EF-9613, EF-9713 contains a 0.035" Nitinol Guidewire





## Catheters for Internal Jugular Approach

### 3 Lumen Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length	Y-Introducer Needle	Quick Insertion Syringe
EG-7623	12Fr x 15 cm curved extensions	60 cm	-	-
*EG-Y7623	12Fr x 15 cm curved extensions	70 cm	-	-
EG-7723	12Fr x 20 cm curved extensions	60 cm	-	-
*EG-Y7723	12Fr x 20 cm curved extensions	70 cm	-	-
EG-I7623	12Fr x 15 cm curved extensions	60 cm	+	-
EG-I7723	12Fr x 20 cm curved extensions	60 cm	+	-
*EG-R7623	12Fr x 15 cm curved extensions	60 cm	-	+
*EG-R7723	12Fr x 20 cm curved extensions	60 cm	-	+

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.035" Marked Guidewire J-Flex Tip
- > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe > 2 Heparin Locks
- (\*) EG-Y7623, EG-Y7723, EG-R7623, EG-R7723 contains a 0.035" Nitinol Guidewire

### 2 Lumen Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length	Y-Introducer Needle	Quick Insertion Syringe
EF-7623	12Fr x 15 cm curved extensions	60 cm	-	-
*EF-Y7623	12Fr x 15 cm curved extensions	70 cm	-	-
EF-7723	12Fr x 20 cm curved extensions	60 cm	-	-
*EF-Y7723	12Fr x 20 cm curved extensions	70 cm	-	-
EF-I7623	12Fr x 15 cm curved extensions	60 cm	+	-
EF-I7723	12Fr x 20 cm curved extensions	60 cm	+	-
*EF-R7623	12Fr x 15 cm curved extensions	60 cm	-	+
*EF-R7723	12Fr x 20 cm curved extensions	60 cm	-	+
*EF-9623	11Fr x 15 cm curved extensions Rigid Tip	60 cm	-	-
*EF-9723	11Fr x 20 cm curved extensions	60 cm	-	-

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.038" Marked Guidewire J-Flex Tip > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe > 2 Heparin Locks
- (\*) EF-Y7623, EF-Y7723, EF-R7623, EF-R7723, EF-9623, EF-9723 contains a 0.035" Nitinol Guidewire

### 2 Lumen Bent Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length	Y-Introducer Needle	Quick Insertion Syringe
EF-7693	11.5Fr x 15 cm	60 cm	-	-
*EF-Y7693	11.5Fr x 15 cm	70 cm	-	-
EF-7793	11.5Fr x 20 cm	60 cm	-	-
*EF-Y7793	11.5Fr x 20 cm	70 cm	-	-
EF-I7693	11.5Fr x 15 cm	60 cm	+	-
EF-I7793	11.5Fr x 20 cm	60 cm	+	-
*EF-R7693	11.5Fr x 15 cm	60 cm	-	+
*EF-R7793	11.5Fr x 20 cm	60 cm	-	+

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.038" Marked Guidewire J-Flex Tip > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe > 2 Heparin Locks
- (\*) EF-Y7693, EF-Y7793, EF-R7693, EF-R7793 contains a 0.035" Nitinol Guidewire

### 2 High Flow Catheter Sets

#### High Flow Femoral Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length
EF-8613	14Fr x 15 cm	60 cm
*EF-F8613	14Fr x 15 cm	60 cm
EF-8713	14Fr x 20 cm	60 cm
*EF-F8713	14Fr x 20 cm	60 cm
EF-8813	14Fr x 24 cm	70 cm
*EF-F8813	14Fr x 24 cm	70 cm

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.038" Marked Guidewire J-Flex Tip > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe > 2 Heparin Locks
- (\*) EF-F8613, EF-F8713, EF-F8813 contains a 0.035" Nitinol Guidewire

### 2 High Flow Jugular Bent Catheter Sets

Art.-Nr.	Catheter Description	Guide Wires Length
EF-8693	14Fr x 15 cm	60 cm
*EF-F8693	14Fr x 15 cm	60 cm
EF-8793	14Fr x 20 cm	60 cm
*EF-F8793	14Fr x 20 cm	60 cm

Qty/Box: 10

Catheter Tray Contents

- > Radiopaque PU Catheter > 0.038" Marked Guidewire J-Flex Tip > 18Ga x 7cm Introducer Needle > 2 Vessel Dilators > 5ml Syringe > 2 Heparin Locks
- (\*) EF-F8693, EF-F8793 contains a 0.035" Nitinol Guidewire