

## For those little ones in need of some peace and quiet...

Designed specifically for premature babies, MiniMuffs noise attenuators protect their sensitive ears, and provide a safe environment for healthy development

Provides a "quiet hour" for premature babies in the NICU

- Helps meet American Academy of Pediatrics guidelines for safe noise levels in the NICU<sup>1</sup>
- Reduces sound levels by at least 7 decibels (dB)<sup>1</sup>
- Results in a reduction of sound pressure level of more than 50%<sup>1</sup>

Ideal during special patient care procedures and transport

- Protects the sensitive ears of babies under ventilation
- Decreases the MRI noise levels when used in conjunction with earplugs<sup>1</sup>
- Reduces disruptive noise associated with transfer of infants to different hospital units and provides quieter environment during transport in emergency vehicles

Safe and comfortable

- Fits comfortably around the baby's ears with a soft foam oval-shaped design and gentle hydrogel adhesive
- Meets infection control standards through single-patient use design

Reduced noise levels are clinically demonstrated to positively impact infant state<sup>1</sup>. By reducing noise levels, MiniMuffs noise attenuators result in potentially longer periods of sleep time, optimal oxygenation and improved physiological stability.



The MiniMuffs neonatal noise attenuators pack includes:

- 36 individually wrapped pairs
- Complete Instructions for Use



### Specifications

#### Nominal Dimensions

<b>Inner width:</b>	1.125 in.	2.856 cm
<b>Inner length:</b>	1.375 in.	3.493 cm
<b>Inner height:</b>	0.31 in.	0.787 cm

**natus**<sup>®</sup>

**Natus Medical Incorporated**  
1501 Industrial Road  
San Carlos, CA 94070 USA  
1-800-303-0306  
+1-650-802-0400  
[www.natus.com](http://www.natus.com)

<sup>1</sup> References on file at Natus

© 2013 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.  
P/N 050312K