

# Nicolet Ambulatory Monitoring



## Nicolet EEG Wireless Amplifier

The Nicolet EEG Wireless Amplifier is ideal for performing ambulatory EEG and home sleep testing procedures giving patients the flexibility to perform their daily activities.

- Allows you to interface with standard Nicolet EEG machine
- No additional software required
- Built in SPO<sub>2</sub> module provides ECG, SPO<sub>2</sub> & Pleth
- 32 channel includes 9 bipolar channels
- DC Input box provides ability to record high level DC ancillary equipment such as EtCO<sub>2</sub> and CPAP
- External batteries permit extended recording
- Batteries are easily hot-swappable by patient
- Compact device and garment design ensures patient comfort

## Nicolet Ambulatory EEG

Easy-to-use Software Applications  
Completely integrated with Nicolet EEG software suite and NicVue patient data management software providing:

- Spike & Seizure
- Mapping
- Sleep
- Trends
- Support Connect
- HIS/HL7 Compatibility

### System features

- Record continuously for up to 72 hours, with data review directly from 2GB flash card, or record in EDF onto an ATA flash disk
- Combined with Nicolet's EEG robust software for review, you can offer a robust service to your patients



Nicolet Ambulatory EEG and Sleep screening effectiveness gives patients the flexibility needed to go about their daily routine while reliably recording data.

[www.natus.com](http://www.natus.com)

## Celebrating 60 Years . . .

