

## The ultimate standard for accurate weight



Olympic Smart Scale Models 60 and 65 combine remarkable accuracy with simple operation. To weigh a baby, simply place the child in the tray, and the scale will automatically compute weight to the nearest gram – even if the baby is moving.

With the Smart Scale, changing numbers do not have to be tracked – the final weight is automatically locked on the display. One can continue caring for the baby and chart the weight when convenient. Even if the baby is removed from the scale, weight remains locked on the display.

- Weigh active babies quickly and accurately
- Four-sided Easy Lift trays are safer for babies while weighing
- Available as table-top or roll-around models



Smart Scale  
Model 65



Smart Scale  
Model 60

## Accurate measurements

- Automatically computes weight to the nearest gram
- Automatically zeros itself before each weighing for accurate reproducible weight measurements
- Zero (tare) button allows extra weight, such as blankets or toys, to be cancelled out to accurately weigh only the baby

## Easy to use

- Fully Automatic – after placing baby on the scale, weight is automatically computed and locked on the display
- Automatic Weight Averaging – when the scale detects motion, it automatically averages and computes within seconds, and locks the final weight on the display
- Memory Recall button – allows display of the last averaged weight, even if the scale has been turned off
- Four-sided Easy Lift Tray (Smart Scale only) – Safer for baby, and hand ports make it easier for lifting baby in and out
- Smooth contours, molded shape and edgeless tray enable easy cleaning

## Efficient

- Weigh babies within seconds, even when they are moving
- Weighs on rechargeable batteries or connected to an electrical outlet while batteries are charging
- Automatic shut down feature conserves battery life
- Available as either a table-top style or as a roll-around model on a sturdy cart

## Roll-Around Smart Scale® Model 65



- Patient is protected with scale securely anchored to the sturdy roll-around base
- Padded crash-bar handles at both ends protects the scale when moved around
- Lower height makes lifting infants in and out of the scale safer for the baby, while reducing potential back strain for the clinician
- Four fully-swiveling, ball-bearing casters enable easy maneuvering
- Supplies can be stored in large enclosed cabinet with transparent doors

## Technical Specifications

### Resolution

1 gram (scale can also be set to read 5 gram resolution)

### Locking Units to Grams (optional)

Scale can be set to measure in grams only.

### Capacity

Infant scales safely hold up to 15,000 grams (33 lbs).

### Battery Operation

Internal Ni-Cad batteries included. Automatically recharges when connected to an electrical outlet.

### Front panel "charging" signal

Scale can be used while batteries are charging.

### Calibration

Easy calibration procedure can be performed by local technician. Standard calibration weight available; see ordering information.

### Electrical Safety

CSA C/US certified to CSA 60601.1 and UL 60601-1.

### Electrical Requirements

120V~, 60 Hz, 0.2A, six-foot power cord with hospital-grade connector.

## Ordering Information

### Olympic Smart Scale Model 60 – Table-Top

Complete with nursery tray and rechargeable batteries.

**PART NO. 56360**

### Olympic Smart Scale Model 65 – Roll-Around

Complete with nursery tray and rechargeable batteries.

**PART NO. 56365**

### Calibration Weight

7500 ±1 gm. Nickel plated steel.

**PART NO. 56315**

**For Information or Ordering**  
**Call 1-800-303-0306**

(toll-free in US/Canada)

**natus**

Natus Medical Incorporated

1501 Industrial Road  
San Carlos, CA 94070 USA  
1-800-303-0306  
+1-650-802-0400  
[www.natus.com](http://www.natus.com)