

RIGID ENDOSCOPES HYSTEROSCOPES

OUTSTANDING QUALITY IN ENDOSCOPIC VISION


- High resolution
- Clear images
- High resistance



OPTOMIC[®]
producing quality

OPTOMIC offers a complete range of rigid hysteroscopes for the inspection of the uterine cavity. They allow for the diagnosis of intrauterine pathologies and serve as a method for surgical intervention. All the OPTOMIC hysteroscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the specialist.


RIGID ENDOSCOPES HYSTEROSCOPES



ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES



PEEK EYEPIECE LENSES

	Ø	Length	0°	12°	30°	70°
HYSTEROSCOPE	2,7mm	270mm	● ●		● ●	
	4mm	300mm	● ●	● ●	● ●	● ●

● Non-Autoclavables ● Autoclavables

- Endoscope with rod lenses system.
- Adapters for connection to type *Storz/Olympus*, *Wolf* and *ACMI* fiber optic cables.
- Body made of stainless steel sealed by laser welding.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.
- Eyepiece of PEEK resistant to heat.

OPTOMIC ESPAÑA S.A.
 Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
 Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
 export@optomic.com • www.optomic.com



producing quality

Specifications are subject to change without notice

