

CHAIR FOR GYNAECOLOGY

OP-G7

OUTSTANDING ADJUSTMENT AND COMFORT



- Compact
- Funcional
- Ergonomic

The OP-G7 is a functional examination and treatment chair for Gynaecology and Obstetrics. It is equipped with all necessary features to meet the current demands for flexibility and ergonomics.

CHAIR FOR GYNAECOLOGY **OP-G7**



FEATURES

UP / DOWN MOVEMENT
 BACKREST MOVEMENT
 LEGREST MOVEMENT
 LEG SUPPORTS MOVEMENT
 TRENDLENBURG
 ELEVATION HEIGHTS
 STAINLESS STEEL BASIN
 PAPER HOLDER
 UPHOLSTERY

FIRE PROOF CLASS
 POWER SUPPLY CONNECTION
 MAXIMUM LOAD
 MEMORIES

OP-G7

Motorized column
Motorized 110° to 190°
Optional motorized
Manual
-10° inclination
595mm min – 1.000mm max
325 x 265 x 65mm
480/490mm 70m roll
PUR-PVC-PO - Washable with
antimicrobial treatment
BS 5852 Part.1 and 2 // M2 (NF P 92503)
110/240v /50-60Hz.
150kg
2 (4 additional optional)



Sturdiness and Stability

Mounted on a telescopic column, which provides a great stability and length range, the OP-G7 chair rests on a quality steel structure which gives it an exceptional sturdiness.

Backrest and seat movement

A linear motor moves the backrest, allowing the specialist to achieve the desired position.

Chair controls

Using the foot control the chair can be moved in every way: up and down and the backrest backwards and forwards. It has up to six direct memories. All these movements can also be activated through a remote control (optional).

Legrest

Electric movement helped by motor turns the gynaecological chair into an obstetrics table.

Leg supports

Completely built in stainless steel, the OP-G7 leg supports can be adjusted to any position, and thanks to their system of self-retaining ball and manual anchorage, they hold the patient's legs firmly. They are padded in flexible polyurethane, very long-lasting and easily washable.

Hygiene

The chair's upholstery material is synthetic, impermeable, washable and specially treated against bacteria (Sanitized), as well as fire-proof class M2 (NF P 92503). Integrated paper roll holder to spread out the cellulose sheet over the backrest and seat. The base's plastic trim can be lifted easily to clean under the chair.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
 Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
 export@optomic.com • www.optomic.com



producing quality



Specifications are subject to change without notice

