

ENT CHAIR

OP-S4

COMFORT AND SIMPLICITY IN THE ENT PRACTICE

- Completely electrical
- Sturdy
- Ergonomic



OPTOMIC[®]
producing quality

The OP-S4 ENT chair has been developed to meet the needs of the modern ENT practice. It combines sturdiness and simplicity in its design. With its ergonomic design, any ENT examination and minor surgical procedure can be carried out in the office.

ENT CHAIR OP-S4

OP-S4

Electrical
Electrical
Foot-operated joystick
12 mm thick treated steel
PUR-PVC-PO - Washable with antimicrobial treatment
Blue standard (choice of other colours)
BS 5852 Part.1 and 2 // M2 (NF P 92503)
Metallic and fold-down
Automatic fold-down with backrest movement
Variable height and inclination. Removable
Adjustable headrest optional
330°
-10° inclination
From 90° to 190°
Max. 900 mm - Min. 550 mm
110/240v/50-60Hz
490mm (w) x 650mm (h)
150 Kg



Fold-down footrest



Control switch



Rotation brake

UP/DOWN MOVEMENT
 BACKREST AND LEGREST MOVEMENT
 MOVEMENT CONTROL
 STRUCTURE
 UPHOLSTERY
 UPHOLSTERY COLOUR
 FIRE PROOF CLASS
 FOOTREST
 ARMRESTS
 HEADREST
 CHAIR ROTATION
 TRENDELBURG POSITION
 BACKREST ANGLE
 SEATING HEIGHT
 MAINS CONNECTION
 DIMENSIONS (BETWEEN BASE TRIMS)
 MAXIMUM LOAD



Ergonomics

The OPTOMIC chair's design is based on a search for the comfort of the patient and ergonomics during examination. This is the reason the geometry of its structure adapts to achieve the optimum distance between patient and doctor in every step of examination and treatment.

Sturdiness and Stability

The base of the OP-S4, although compact, gives it great stability and sturdiness, making it appropriate even for large patients.

Strong Armrests

They are composed of a stainless steel interior, and padded with a flexible and washable material. They fold down automatically, synchronized with the backrest movement. They can also be folded up to get out of the way when the chair is in a sitting position.

Footrest

Fold down and robust. The padded side becomes part of a comfortable and broad examination table.

Backrest and Footrest Movement

The backrest and footrest movements are mechanically synchronized, so any position from 90° to -10° (Trendelenburg) can be achieved using the pedal.

A Precise System

The system is controlled through a foot switch which permits easy handling of the chair. It operates both the up/down movement and the inclination of the backrest.

Rotation

The OP-S4 chair allows for a rotation of 330°.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
 Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
 export@optomic.com • www.optomic.com



Specifications are subject to change without prior notice

