

# HD INTEGRATED ENDOSCOPY SYSTEM STROBOLUX III LED+HD

- HD endoscopy camera
- Stroboscope
- LED light source



The Strobolux III LED + HD integrated endoscopy system is revolutionary in the ENT practice. It allows the user to carry out stroboscopy and continuous light examination synchronized with the Optomic HD camera.

# INTEGRATED SYSTEM STROBOLUX III LED+HD


## HD ENDOSCOPY CAMERA

- HD output, HDMI type.
- 1080i HD digital resolution.
- 2 outputs: S-Video and composite video.
- Shutter speed 1/100,000 sec.
- Manual and automatic white balance.
- Automatic synchronized ELC.
- Synchronization of camera and stroboscope to optimize the image quality.





## STROBOSCOPE

- LED stroboscopic system.
- 2 independent and simultaneous light outputs.
- Phase and slow movement function.
- Analysis frequency from 60 to 1,000Hz.
- Phase range variation from 0° to 400°.
- Voice analysis: spectrogram and waveform.
- Air and contact microphone.
- 5" TFT display.



## LED COLD LIGHT SOURCE

- LED type illumination.
- 2 independent outputs.
- Alternating or simultaneous function of the outputs.
- Life of LED lamp: 60,000h.
- Intensity of use memory.
- Adaptors for all cable brands in the market.



Specifications are subject to change without notice

